

TNFSF12-tumor necrosis factor (ligand) superfamily, member 12 Gene

Reference: PTXBC019047

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APO3L; DR3LG; TNLG4A; tumor necrosis factor ligand superfamily member 12; APO3 ligand; APO3/DR3 ligand; TNF-related WEAK inducer of apoptosis; tumor necrosis factor (ligand) superfamily, member 12; tumor necrosis factor ligand 4A; tumor necrosis factor su		

For research use only.

GENE DESCRIPTION	tumor necrosis factor (ligand) superfamily, member 12
GENE ID	8742
ORIGIN SPECIES	Human
ACCESSIONS	BC019047
VECTOR	pDONR223
NCBI SYMBOL	TNFSF12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgcccgtcggagccagagggcggagggggcgggggggagccgggcaccgcccgtctggtcccgctcgcgtggcctgggctggcgtggcctgctcctggcctcctgctggtccgtgttcagtttgggagccgggcatcgtgtccgccaggagcctgccagga
ggagctggtggcagaggagaccaggacccgtcggaactgaatccccagacagagaagaagccaggatcctgcgcctttcctgaaccgactagtctggcctgcagaagtgcacctaaggccggaaaaacacgggctcgaagagcgtcgcagccattatgaagtcatccacgac
ctggacaggacggagcgcaggcagatggaggttacacaactgtctgaggccatga