

## SP110-SP110 nuclear body protein Gene

Reference: PTXBC019059

|                |                                                                                                                                                                                                                                                             |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | transcriptional coactivator Sp110; sp110 nuclear body protein; IFI41; IFI75; IPR1; VODI; interferon-induced protein 41, 30kD; interferon-induced protein 41/75; interferon-induced protein 75, 52kD; phosphoprotein 41; phosphoprotein 75; speckled 110 kDa |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | SP110 nuclear body protein                                                                       |
| <b>GENE ID</b>          | 3431                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC019059                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SP110                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgttcacatgacagagccatggaagagctcttttcagcactcatgaccagaagctggggatgcctatgccatacacaagccattccctctttgaaggcccttagacaactccatcatcactaagagaatgtacatggaatctctggaagcctgtagaatttgatccctgtatcca  
gagtggtgcacaacattctcacccactggagagagcttttaacctgtctctctgtgacattgttcagtcacaatttaacctgcgtgaatatcccaatctggtgacgatttacagaagcttcaaacgtgttggtctctctatgaacggcagagcagagacacccaatcctactgaagccccc  
ctggcctagcagaaggaagctccctccatacccaactggcgtgccccaccacaaccccccaaccaagctgttcacccgtgtgcgccaagagtcagtgagcctggaacatcctccagcaaacgcatgagatcctgagtgagtcgcccagcccatctgacctgtcctgtcctctccctg  
cactatccaggaaggaagcacttcagtgaccaatgacaagttaacatcccaaatgaatcggaagaagactcagaagagatgccagcctcctcactagcactgtgcaagtggccagtgacaacctgaccccccaataagagataaagaagccctcaagagatgccccactct  
cccttggctctatgccagagataagagataattccagaacaaatgaccagaagagcccgagggtgtccagcacaccttcagacaagaaggaagaaaagaaaaagatgtatctgtgcaactccaaaaaggagacataaagaaaaagcctcccaagagccttatctagaca  
cggaatccaaaagaagctcaaaagggtggtcaggttctcaaaaagaagatgactcaactgttaactccacgttagagacaaagggcccaaaaggcgagaactgaatgtgcccgaagtcgagatcagagagatcattgatggcacttcagaatgaatgaaggaagaggtccag  
aagacgctagtacaccacgaagggtcacacaagggcgacgtcactgggcatggcactcaagaagaagctccaagtggtggaaggtgactcaaggaagacgactcaacctggaactcagaggtcatgatgagggtccaaaaggcaagaactaatgtgccccgaagtcag  
atcgaaagaaaaaggaagagatctgttcaagctcaaaaagggaatttcagaaaaatattaccgaagaggaagaaaccccaaaagtgacactgtggaatttctactgttcaagctccccgtgacctgtgtgtgagcgaaagggaatttatataagaagaaatgaacacggatcct  
cagtgaaagtcattcggaatgagatggaactgtgtaacaccaaataatgaattgaagtcgaaggaaaagggaaggaacgcaagaactggaacggaatatactgtggaaggaaatgacccctaggagagctgctgaagagtggaactttgtctgtcctccaaagataatctcaagagaga  
gttaataagcaagtga