



Reference: PTXBC019247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BS-84; CILD28; CT142; HSD-1; SMP1; SPAG3; hSMP-1; sperm-associated antigen 8; sperm membrane protein 1; sperm membrane protein BS-84; testicular tissue protein Li 177; sperm associated antigen 8		

For research use only.

GENE DESCRIPTION	sperm associated antigen 8
GENE ID	26206
ORIGIN SPECIES	Human
ACCESSIONS	BC019247
VECTOR	pDONR223
NCBI SYMBOL	SPAG8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr