



Reference: PTXBC019248

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	monoacylglycerol lipase ABHD2; HS1-2; LABH2; PHPS1-2; 2-arachidonoylglycerol hydrolase; abhydrolase domain-containing protein 2; alpha/beta hydrolase domain containing protein 2; lung alpha/beta hydrolase 2; protein PHPS1-2; testicular tissue protein Li 6		

For research use only.

GENE DESCRIPTION	abhydrolase domain containing 2
GENE ID	11057
ORIGIN SPECIES	Human
ACCESSIONS	BC019248
VECTOR	pDONR223
NCBI SYMBOL	ABHD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]