



Reference: PTXBC019257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMT2Z; ZCW3; ZCWCC1; MORC family CW-type zinc finger protein 2; zinc finger CW-type coiled-coil domain protein 1; MORC family CW-type zinc finger 2		

For research use only.

GENE DESCRIPTION	MORC family CW-type zinc finger 2
GENE ID	22880
ORIGIN SPECIES	Human
ACCESSIONS	BC019257
VECTOR	pDONR223
NCBI SYMBOL	MORC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]