

NHP2L1-NHP2 non-histone chromosome protein 2-like 1 (*S. cerevisiae*) Gene

Reference: PTXBC019282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NHP2L1; 15.5K; FA-1; FA1; NHPX; OTK27; SNRNP15-5; SPAG12; SSFA1; NHP2-like protein 1; NHP2 non-histone chromosome protein 2-like 1; U4/U6.U5 tri-snRNP 15.5 kDa protein; [U4/U6.U5] tri-snRNP 15.5 kD RNA binding protein; high mobility group-like nuclear pro		

For research use only.

GENE DESCRIPTION	NHP2 non-histone chromosome protein 2-like 1 (<i>S. cerevisiae</i>)
GENE ID	4809
ORIGIN SPECIES	Human
ACCESSIONS	BC019282
VECTOR	pDONR223
NCBI SYMBOL	NHP2L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgaggctgatgtgaatccaaaggcctatccccctggcgatgcccaacctcaccagaagctactggacctcgttcagcagtcgtgaactataagcagcttcggaaggagccaatgaggccacaaaacctcaacaggggcatctctgagttcatcgtgatggctgcagacgccg
agccactggagatcattctgcacctgccgctgctgtgaagacaagaatgtgccctacgtgtttgtgcgtccaagcaggccctggggagagcctgtgggtctccaggcctgtcatcgctgttctgtccatcaaagaaggctcgcagctgaacagcagatccaatcattcagcag
tccattgaaggctcttagtctaa