



Reference: PTXBC019337

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin heavy constant gamma 1 (G1m marker)
GENE ID	3500
ORIGIN SPECIES	Human
ACCESSIONS	BC019337
VECTOR	pDONR223
NCBI SYMBOL	IGHG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]