

TSPAN18-tetraspanin 18 Gene

Reference: PTXBC019342

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSPAN; tetraspanin-18; tspan-18; tetraspanin 18		
For research use only.			

GENE DESCRIPTION	tetraspanin 18
GENE ID	90139
ORIGIN SPECIES	Human
ACCESSIONS	BC019342
VECTOR	pDONR223
NCBI SYMBOL	TSPAN18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggcgaactgtctgagctcatgaagtatctgatgtttgtattcaattttctatatttctggcgggggcctgctgctggccatcggcatctgggtcatgttggaacccaccggcttcggggagatcgtggctgccaatcctctgctctcacggcgccatcctctggccatggg
gggctgtcttttctgctcggcttctgctggctgctgctggggcctgctgagacaagtgtctgctgtatttttctctgtttcatcctgatcatcttctggcagagctctcagcagccatcctggccttcatcttcaaggaaatctcaccgagaattttcaccaggagctcaccagcact
accaggccaataacgacacagacttctctgcccactgggaactcgggtcatgatcatttgggtgctgctggggcctcaacgggctgaagactttaagtttgcactgtgttttgcactcctgacctggatagtgagagggtgcccggggcctgctgcccggagggaaccccaagtgcgga
cggggtcctgctgagccgggaggagtgctcctcctgggaaggagcctatttctaaacaagcagggtgttacacgggtgatcctcaacacctcagagactacgtctacttggccggagccctgcccactggggtactggccatcgaggagaagaagagctcattctgtctcacaagccaccgg
catcctgtatcagcttcagcctccctcaggcttccagtcaccagggaactcggagccacagcctag