

TMEM174-transmembrane protein 174 Gene

Reference: PTXBC019346

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 174		
For research use only.			

GENE DESCRIPTION	transmembrane protein 174
GENE ID	134288
ORIGIN SPECIES	Human
ACCESSIONS	BC019346
VECTOR	pDONR223
NCBI SYMBOL	TMEM174
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcaggcagcggccgttgaggactccctgtcaatgtgttcctcactccttacacaccagcaccgctgacatccaggtgccgatgatgacaaggcggggccacctgtctcttcaggcatctttgggaactgtgtgggacacattcactgtcatgggtggatcaaa
taccagggtgtctccacttgaatggaccagctccttgggcccgtcctgtcagttgggtgacattcctctgattgtgtgtgcaagttcaaaatgctcctgccagttgtgcaagaaagtggagaaagggtcccggactcggacagacaccaggaggaccatcattgtttcac
tggcatcaaccaaccatcacttccatggggccactgtgtgtagtacctcctccttatgggttcacagcctatggggataaataccagctacctgcagctgtgtgtgagccctgcggcctataacctctggagggcagcagcccatgtcaagtctcctcaatactacac
catctaccctcaagataactctgcatctgtgtgtgatgagggtgcctttcttcacggacgttgaaatcacaggcccaatcctgatgttgaccagctagaagagacacagctggagagaggcctgtgcctgtctctcctcccttatgaagaataatactctcctcctcctag