

FABP5-fatty acid binding protein 5 (psoriasis-associated) Gene

Reference: PTXBC019385

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E-FABP; EFABP; KFABP; PA-FABP; PAFABP; fatty acid-binding protein, epidermal; epidermal-type fatty acid-binding protein; fatty acid binding protein 5 (psoriasis-associated); psoriasis-associated fatty acid-binding protein homolog; fatty acid binding prote		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 5 (psoriasis-associated)
GENE ID	2171
ORIGIN SPECIES	Human
ACCESSIONS	BC019385
VECTOR	pDONR223
NCBI SYMBOL	FABP5
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccacagttcagcagctggaaggagatggcgcctggggacagcaaaagccttgatgaatacatgaaggagctaggagtggaatagctttgcgaaaaatggcgcaatggccaagccagattgatcatcactgtgatggtaaaaacctcaccataaaactgagagcacttga
aaacaacacagtttctgtacctgggagagaaagttgaagaaccacagctgatggcagaaaaactcagactgtctgcaactttacagatggtgcattggttcagcatcaggagtgatgggaaagaaagcacaataacaagaaatgaaagatggaaatgagtgagtggtc
atgaacaatgtcacctgtactcggatctatgaaaaagtagaataa