

C18orf51-chromosome 18 open reading frame 51 Gene

Reference: PTXBC019628

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C18orf51; protein FAM69C; family with sequence similarity 69 member C		

For research use only.

GENE DESCRIPTION	chromosome 18 open reading frame 51
GENE ID	125704
ORIGIN SPECIES	Human
ACCESSIONS	BC019628
VECTOR	pDONR223
NCBI SYMBOL	C18orf51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctttttgaacctaaatgagggaaatccttgagcaaaactgcacaggagatgaagactgcaatttcttgactgttttcaagatgtgattacgagtcacaaatgcggagcgcagcgcgtaaacaacaacctgcaggtcatctgtgacaaaatattcgccattggtttccgcgcctct
caagagctctgcggctctttccagcttcagctgcagttacaggagggcgtgcaggaaatgtgcagacctgggggtccccagtggaacacccggagagcagcctccagcgtgttctggaagcttcgccaactcctccaagccacactgaggagctgcaggaggcagagaagtag