

## SDR42E1-short chain dehydrogenase/reductase family 42E, member 1 Gene

Reference: PTXBC019670

|                |                                                                                                                            |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | HSPC105; short-chain dehydrogenase/reductase family 42E member 1; short chain dehydrogenase/reductase family 42E, member 1 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | short chain dehydrogenase/reductase family 42E, member 1                                          |
| <b>GENE ID</b>          | 93517                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC019670                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SDR42E1                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggaccacaaagatctcaaaaggaaagtgtcctcattacaggaggagtggtctatgttgcctgggctgtgccctgaacaaatggagtcattgttggacatcagcagccctgtcctcaaccattcagaaggaatcaagttatacaaggagacatccgccctgtc  
tgacgtagagaaagcctccaggtatgcagacgtcactgtgtgtccatattgctcttattggtatgtcaggggcgggagcaactcaatcgaacctgatcaaaagatcaacgtcaggggcacagacaacatcctcagggttgcctcaaaaggagaagggtgccaggttagttacaccagca  
cttcaatgtcatcttggaggtcaagttatcagaatggggtatgaatctctgcccctacctgctctcactccacccctgatcactactctcggacaaagtcaattgcagagcagaaggtgtcggaggcgaatgctacacccctggacagaggcgacggtgtcttaagaacctgcgctctgag  
gccagctggcatctatgggctgggagacaagacacctccaggatagtcagctacatcgagaagggtctgttcaagttgtctacggggacccaggagcctggttgagttgtccacgtggataacttggtgcaggtcactctcctcagagccctgagagctgacaagggc  
catattgctctgggagccctactctcatctcagatggcagaccgtgaacaactttgattctccggcctctggttgaggccctgggctacacattcccgtctaccggcctgacattgacattggtctactgcttgccttttaacagagatggtcactcatttgggtcgactctacaactcc  
agcccttctcactcgactgaagtttacaactgggtgtcacacattatttagcttagagaaagccaagaagagctaggttataaggctcagccattgacctccaggagcagtggaatggtttaaagccatggctatggcagaagttctggaagtcgtgactcggaggtgtttgttggg  
atgggctattgtcttctcctcctgattatagcagttctcatgtg