

C12orf49-chromosome 12 open reading frame 49 Gene

Reference: PTXBC019843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0454 protein C12orf49 homolog; c12orf49; chromosome 12 open reading frame 49 S homeolog		

For research use only.

GENE DESCRIPTION	chromosome 12 open reading frame 49
GENE ID	79794
ORIGIN SPECIES	Human
ACCESSIONS	BC019843
VECTOR	pDONR223
NCBI SYMBOL	C12orf49
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaacctggcgccatggtgtggcgccgcttctgcggaagaggtgggtgctcgccctggtctcgggctgtcgtctgtctactctcagcagcaccttcaagcaggaggagggcagtgagagataggaatctctcaggttcatgaccataatcagcccatccgtggaaa
gtgcagtttaacttgggcaatagcagtcgtccgagcaatcagtgccgcaactcattcaagggaagcacctcatcacggatgaactcggtctgttgcgagagggaagatttgcgtgtaaatggctgtgtaatgtcaacgtccctagcaggaagcagtgctgtgatggctgtgccc
caacggctgtgcagcgctatgagtactgtgtctctcgtcgtgcagcccaacaagcaacttctcctggagcgttctcctcaaccgggcagccgtggcattccagaaccttctatggcagtcgaagatcattgagttgtgcctggccaaatgcaggacctatctcagagcgtgcagc
atgagaacacctaccgggaccccatagcaaatgattgctatggagaaagcccgccgagcttctcccgcttga