

TMEM42-transmembrane protein 42 Gene

Reference: PTXBC019851

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 42		
For research use only.			

GENE DESCRIPTION	transmembrane protein 42
GENE ID	131616
ORIGIN SPECIES	Human
ACCESSIONS	BC019851
VECTOR	pDONR223
NCBI SYMBOL	TMEM42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgagaggccgggctccggcgccgtgtccgacccgctgacaccccgcggaattccctccgcacctccaggcgggtgcgatgcggcgccgttttggggcgattcaactgtgtgtgcgcggcggttggggccctggccgcccctccgccaagct
ggccttcggcagcgaggtgagcatgggtttatgcgttaggcattattgtgatggcgagcacaattctctgatgtggaccttcttagccggggcctcagtttctccatgtctcagccattgcattctgtcacagtgcattttcaaatatcctcagctcggccttctgggctatgtgtgtatgga
gagtgccaggaggctctgtgtggggaggaggtgttcttattctctgcgactcaccctaataccacagggaagctccacccacctggaagcccttccacacaagcagcagtag