

CLEC2D-C-type lectin domain family 2, member D Gene

Reference: PTXBC019883

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLAX; LLT1; OCIL; C-type lectin domain family 2 member D; C-type lectin related f; C-type lectin superfamily 2, member D; LLT-1; lectin-like NK cell receptor; lectin-like transcript 1; osteoclast inhibitory lectin		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 2, member D
GENE ID	29121
ORIGIN SPECIES	Human
ACCESSIONS	BC019883
VECTOR	pDONR223
NCBI SYMBOL	CLEC2D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatgacagtaacaatgtggagaaagacattacaccatctgaattgcctgcaaacccagggtgtgtgcatcacaagagcattctattaaagctacctaatttggcgcttattttcttaacatgtttctgacaatcatagtgtgtggaatggttgcctttaagtcaataagagctaactgcca
tcaagagccatcagatgtcttcaagctgcatgccagaaagctggattggttttcaagaaagtgtttctattttctgatgacaccaagaactggacatcaagtcagaggtttgtgactcacaagatgctgattctcagggtgaaagctccaggaaactgaatttcctgttgagatataaagg
ccatctgatcactggattggcctgagcagagaacaaggccaacatggaaatggataatggtactgaatggacaagacagtaa