



Reference: PTXBC019944

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLVR1; Glvr-1; PIT1; PiT-1; sodium-dependent phosphate transporter 1; gibbon ape leukemia virus receptor 1; leukemia virus receptor 1 homolog; solute carrier family 20 (phosphate transporter), member 1; solute carrier family 20 member 1		

For research use only.

GENE DESCRIPTION	solute carrier family 20 (phosphate transporter), member 1
GENE ID	6574
ORIGIN SPECIES	Human
ACCESSIONS	BC019944
VECTOR	pDONR223
NCBI SYMBOL	SLC20A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]