

## TBCC-tubulin folding cofactor C Gene

Reference: PTXBC020170

|                               |                                                               |                     |                |
|-------------------------------|---------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                           | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                          | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | CFC; tubulin-specific chaperone C; tubulin folding cofactor C |                     |                |
| <b>For research use only.</b> |                                                               |                     |                |

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tubulin folding cofactor C                                                                        |
| <b>GENE ID</b>          | 6903                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC020170                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TBCC                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggagtcctcagttgctccgctgctgctcaggaccggagacatggagtcgccagcgggacctgagcctggtgcctgagcggcttcagagacgcgaacaagcggcagctggaagtgaaggcggaaacaaagcggcagaaccaggaggtagagaaggagaacagccac  
ttttcgtcgcacctttgctcgggagcggagcggcgtggaagagcttctggagcgcgcggagtcggtcggcggctggaaggagcggcctcgcctccaggggctgcagaaactaatcaacgactcagtttttctagccgcttacgacctgcggcagggacaagaggcgtggc  
gcggctgcagcgcggccttggccgagcggcggcggcggcgtgcagcccaagaagcgtttcgtttcaagaccggggaagagatgctgcttctaccaaagtagacgcggctcctggcatccccccggcagttgaaagcatacaggactccccgtgcccaagaaggcggaaagga  
gacctggcccccagctgggtctgctggttttccaactggagtcaccaagcttggagaagagaccagcaggtgacaccagcgcgacgtttttgaccgaactgagcaactgcacggtcagactgtatggaatcccaacacctgcgggttaaccaaggccacagctgcaagctgctc  
tgcggtccggtgctacctctgtttcctggaggactgcagtgactgctgctggcagtgccctgccaacagctcgcatacagtaacgaaagacaccgcacttctcctgaggtgaccagcaggccatcgtggaggactgcagtgaggatccagttcgcccttacactggagctacc  
cggagatcgacaaggacttcagagctcgtgttagataggagcaaaaataactggaacgatgttgacgattttaactggctggccgggatatggcctcccaactggagttattctcctgaaggagcgaataatccagtgaggactaa