

CALN1-calneuron 1 Gene

Reference: PTXBC020200

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CABP8; calcium-binding protein 8; calcium-binding protein CABP8; calneuron I		
For research use only.			

GENE DESCRIPTION	calneuron 1
GENE ID	83698
ORIGIN SPECIES	Human
ACCESSIONS	BC020200
VECTOR	pDONR223
NCBI SYMBOL	CALN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgttccacatgtgaccgccgcttgtgtacaaggggaattacccaacgacgctctctgctgagcagcgaacagctggcctaataatctcctggaggagctcgatgaaatccgagaggcctttcgggttctggaccgggatgggaacggcttcattccaagcaggagc
tgggcatggccatgcgctctttgggtacatgccaagcgaggtggagctggccatcatcatgcagcgcctggacatggacgggatggccagggtgattttgatgaattcatgaccattctggcccaaacctgggtcttcagaaggctcgcgatggtttcttgggaacacgatagacagca
tattctggcagtttgacatgcaaaagataactctggaagagttgaagcacattctctatcatgcctccgagaccacctaacgatgaaggacattgagaacatcattatcaatgagggaagagagcctgaatgagacctcgggaactgccaacagagttgaggagtgcatccagaagc
agaacagacagacctgcgtccggaagagcctcatalgcgcctttgctatggccttcacatcagtgctcatgctgattgcagccaaccagatactccggagcgcatggagtag