

MRPL2-mitochondrial ribosomal protein L2 Gene

Reference: PTXBC020212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-22; MRP-L14; RPML14; 39S ribosomal protein L2, mitochondrial; L2mt; MRP-L2; mitochondrial ribosomal protein L2		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L2
GENE ID	51069
ORIGIN SPECIES	Human
ACCESSIONS	BC020212
VECTOR	pDONR223
NCBI SYMBOL	MRPL2
LENGTH	Proteogenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggccctgtgcgactgacccgcgctctgctctgaacctgggccccaccgctcgcgccccccccgagctgttccccgcccccagatgatgaacaatggcctcctccaacagccctgcttgccttgccttcttactctgtgcccctaatgcc
aaccttgtctctggaagctgtaaccaagfacaccattacaccagtgaaatgaggaaagctggggccgagaccacacagccgaatccgggtgcatggtattggcggggccacaagcaacgttatgaaatgacttctcgtttccgctgaggagaccaagtcaggacc
tttgaggagaaggtatccaagtccgctatgatccctgtaggtcagcagacatgctctggttctgctggggcagccggaacgctggatcgcacagaaaacatgcaggctggagatacaatctgaacttaaccacatagccgaatggcagttgctcgggaaaggatgcg
catcctctgggctctgctgtgggacctcatcaaacgtggaagtggagccgggggtgccaatatccgagctcagggagctgtgggtgctactgcggaaggtgaatggcacagccattatccagctccctctaagaggcagatcaggtgctggaacgtgct
agcaacagtagccgagatccaactgatacaaacgggtcattggcaaggcaggtcgaaccgctggctggcaagagcctaactggcgggtgaccgcaagggggctggctggccgaaagattggccactacccccataagagttactgaaagctgctt
ctgcttccccaaagctga
```