



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LOC283693-actin, gamma pseudogene Gene

Reference: PTXBC020234

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

For research use only.

GENE DESCRIPTION	actin, gamma pseudogene
GENE ID	283693
ORIGIN SPECIES	Human
ACCESSIONS	BC020234
VECTOR	pDONR223
NCBI SYMBOL	LOC283693
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattcagatcatgtttggggccttcacaccccgatcaagtacaggcgccattcaggctaccctgtccctcatgcctctcgtcagtgccattggcaccattgggggactcaggcgatggggtcaccacactgcccatcccaagatagcaccctctcccatcatcctgtttggcc
ctggctgcccgggacctggccactactgtgtgaagatctctcgaagggccgctacagcttcatcaccagatgaagtgggggattgtgcgtgacatccaggaaaagtatcgtgtgttccttggaacttcagcagagagatgccatgccgtgtctcttctctctcggagaagagcta
caagctctccacagacctgtgtgtgctccatcatcaagatgagtggttctgtgtccagaggcgctgtccagccctctcccgatggaaactctgtcattccaggagagaccatctgacctccatcatgaagtgtgatgtggatatctcaggggacacatgacacaaatgtgtgtctgtcgc
acattgctcctcagcccaagatcatcagctgatg