

IGLV4-3-immunoglobulin lambda variable 4-3 Gene

Reference: PTXBC020236

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	immunoglobulin lambda variable 4-3
GENE ID	28786
ORIGIN SPECIES	Human
ACCESSIONS	BC020236
VECTOR	pDONR223
NCBI SYMBOL	IGLV4-3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgggtctctctactactgccctcaatttctccacaggtctctgtgtctgtctgtctgactcagccccgtctgcatctgcttctctgggagcctcgaagctcacctgcacccaaagcagagagcacagctcctacacatcgaatgggtatcaacagagaccaggagggtc
ccccagtatataatgaaggttaagagtgatggcagccacaacaaggaggacgggatccccgatcgttcatgggtccagttctgggctgaccgtacctcacttgcacacccagctgacgatgagctgaatatcactgtggagagaccacgattgatggccaagtcggt
tgggtgtcggcgaggaggaacagctgaccgtcctaagtcagcccaaggctgccccctcggtcactctgttccgccccctctctgaggagcttcaagccaacaaggccacactgggtgtctcataagtgacttctacccgggagccgtgacagtgccctggaaggcagatagcagccc
cgtaaggcgggagtgagaccacacccctccaacaagaacaacaagtagcgggccagcagctacctgagcctgagcagtggaagtccacagaagctacagctgcaggtcacgcatgaaggagcaccgtggagaagacagtgggccctacagaatgttca
tag