

CYB5D2-cytochrome b5 domain containing 2 Gene

Reference: PTXBC020263

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	neuferricin; cytochrome b5 domain-containing protein 2; cytochrome b5 domain containing 2		

For research use only.

GENE DESCRIPTION	cytochrome b5 domain containing 2
GENE ID	124936
ORIGIN SPECIES	Human
ACCESSIONS	BC020263
VECTOR	pDONR223
NCBI SYMBOL	CYB5D2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgttgaggtcggagccgtggcctttgtggcctggtgtagccgcagcagcggtaatggcagcacgcttatggcgtgtgggtcccgcgctggccttccataccgaggagctgtctcgtaccgcggccaggggaccggcctgtacttggcgttgc
ggccgtgtctacgatgtgtctccggccggaggaactacgagcctgggtcccactatagcggcttcgagccgagagacatccagagcttctgaccgggactgttctgaagcaggcctcgtggtgatgacgtatccgacctgtcagccgtgagatgctgacactcacaattggctt
cattctatgagaagaattatgtgtgttgggaggtgacagcaggttctacggagagatggcgtcccacccggcactgaccaggtagaagctgcatcaccagagcttgaggccaacaactacagctgcaagagaagcagacattccgcctgcaacggagtgga
gctcagccagggcagccgctcgtgctccagaagatgaggtgtgagcagagactggatggcgtcccaggaaagctgataagccaggtgctaaggagcccgctgctgtgtgtgagaaccaccggccccctagtggccagatcgggacaacctcacagaaat
cgtgggacctggaccaccaacttggcagagtacagagctgccaccgctagccatcacatgctcttccactctaa
```