

RPS14-ribosomal protein S14 Gene

Reference: PTXBC020515

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EMTB; S14; 40S ribosomal protein S14; emetine resistance; ribosomal protein S14		
For research use only.			

GENE DESCRIPTION	ribosomal protein S14
GENE ID	6208
ORIGIN SPECIES	Human
ACCESSIONS	BC020515
VECTOR	pDONR223
NCBI SYMBOL	RPS14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacctcgaagggaaggaaagaaggaagaacaggtcatcagcctcggaacctcaggtggctgaaggagagaatgtatttggtgtctgccatatttgcacacctcaatgacactttgtccatgtcactgatctttctggcaagaaccatctgccgtgtgactgggtggatgaagg
taaggcagaccgagatgaatcctcaccatgtctgtgtgtggctgccaggtgtggccagaggtgcaaggagctgggtatcaccgccctacacatcaaacctccgggccacaggaggaataggaccaagaccctggacctggggcccgatcggcctcagagccctggcc
gctcgggtatgaagatcggcggtattgaggatgtcaccatccctctgacagcactcgcaggaagggggctgcctgtggtcgcctgtctgta