



Reference: PTXBC020565

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	24p3R; BOCT; BOIT; NGALR; NGALR2; NGALR3; hBOIT; solute carrier family 22 member 17; 24p3 receptor; NGAL receptor; brain-type organic cation transporter; lipocalin-2 receptor; neutrophil gelatinase-associated lipocalin receptor; potent brain type organic		

For research use only.

GENE DESCRIPTION	solute carrier family 22, member 17
GENE ID	51310
ORIGIN SPECIES	Human
ACCESSIONS	BC020565
VECTOR	pDONR223
NCBI SYMBOL	SLC22A17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr