

SLC39A3-solute carrier family 39 (zinc transporter), member 3 Gene

Reference: PTXBC020571

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZIP-3; ZIP3; zinc transporter ZIP3; solute carrier family 39 (zinc transporter), member 3; zrt- and Irt-like protein 3; solute carrier family 39 member 3		

For research use only.

GENE DESCRIPTION	solute carrier family 39 (zinc transporter), member 3
GENE ID	29985
ORIGIN SPECIES	Human
ACCESSIONS	BC020571
VECTOR	pDONR223
NCBI SYMBOL	SLC39A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaattgctagtgccaaatcctgtgcatgggtggcggtgtttctttcatgctgctggctccctgctccccgtgaagatcatcgagacagattttgagaaggcccatcgctcgaaaaagatcctctctctgcaacacctttggaggagggtgtttctggccacgtgttcaacgctctgctccccgtgtgagggaaggtaagggtcctctggcactagcagcagccctggcaccttatggccaaggactctgatgcatttcaacactgatgccaagtgcagtgaaggccttgatgctgtaa