

DDX60-DEAD (Asp-Glu-Ala-Asp) box polypeptide 60 Gene

Reference: PTXBC020601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEAD (Asp-Glu-Ala-Asp) box polypeptide 60; DEAD box protein 60; DEAD-box helicase 60; DExD/H-box helicase 60		

For research use only.

GENE DESCRIPTION	DEAD (Asp-Glu-Ala-Asp) box polypeptide 60
GENE ID	55601
ORIGIN SPECIES	Human
ACCESSIONS	BC020601
VECTOR	pDONR223
NCBI SYMBOL	DDX60
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaatattggaggactttaccactttctacgaattgtttccaaactggctgatatgaatcaggaatatcaactccattgtcaaaaatcaaatcacaggtaaagaatgtgaagactctcaactcgtatctcatttgatgagctgaaggaaaggagtagcaatttcaccattgtttgtctgt
ctgggaactttgatgatgtttgcttcgactagaactccaaccatgttactctagtcacaatcggtgtcaatcgtctcaggtccagtgtgtgtcacagaaattgataaccgaggaaggaaaatgtcgttaatgcctatgcactggatttctacaacatgggttccttgataggattagtc
aggataacaggatgaatgaaggagatgcttattttgtgaaggattttgcactcaccattaaatctatcagtggttccttgcgtgagctatgtgaaaatgaagacgacaacgttgcttagcctttgaacaactgagtacaactttttggaaaagttaacaaagtctaa