

C10orf97-chromosome 10 open reading frame 97 Gene

Reference: PTXBC020605

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf97; CARP; DERP5; MST126; MSTP126; my042; ubiquitin carboxyl-terminal hydrolase FAM188A; CARD-containing protein; caspase recruitment domain containing pro-apoptotic protein; dermal papilla-derived protein 5; protein FAM188A; family with sequence sim		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 97
GENE ID	80013
ORIGIN SPECIES	Human
ACCESSIONS	BC020605
VECTOR	pDONR223
NCBI SYMBOL	C10orf97
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgaactgactaaagagctgatggagctggtgtgtgggaccaaagagcagccccgggtctctcggacacacatttctgccgctggacgcaaggggtgtgttagtgatcagagggatctgcattagaacagttgaaggtggccccgtgtgttattgcacgtgttcaggcattctttg
aagaagctcctgtttcttcggagaagctctcttggcgggattgttcagaggaagagcagaaggaactccttctgacacattgtgtgatatttagaaggtgctgtgttgaccactctggatcactactgctgtttcatgtgttaagaggaaagacaactgaggaaactgctagtattctgggagt
cctgcagagctagttgccaagtgggaacattctaa