

CXorf21-chromosome X open reading frame 21 Gene

Reference: PTXBC020611

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein CXorf21; 5430427O19Rik; chromosome X open reading frame 21		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 21
GENE ID	80231
ORIGIN SPECIES	Human
ACCESSIONS	BC020611
VECTOR	pDONR223
NCBI SYMBOL	CXorf21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcagaagggtatctcagtgactgagtggaatgacatccatggagttgtgcctttataatgagcaggtggctgggaaaggaagagagacaaattctgttgcacccttctattctctgttgatgaaacaaagtcagaagtcctacgtgagctgcaaatcatctgg
caagtttatctctcagtgattcaagagagagccaacatagcagaagtcagagagtcacagtgctgcagacaaacccaatcctgtgtttgaaagcccaacttggctgcagttgaaatagtagagatgccagcagagagacactacttggttccatctcttgcaaaagtattgcaagaatt
ataatgacttacagattgcagggggccaggtgatggccattaatcagtgacacagattttccctctgagagcagttttgaatagggccctttgtcgaagtcatctgagattcctttacccatggaggattccatttctacagccagtgactttcccaaaacatccagcgggtactatcct
attggagaatacaagcatcaaaagaaaaagcagcttcgaatgcagaatcctatttctaatgcagtgctgaatgagtagctggagcagaaggtgttgaggttatataaaacagtacattatggacaccgtgtttcatgacagttctctaccagattctgcgtctgaactcatcatgacaagtg
tagaccaatcagtcctcaagtgctagagagaagaatctggagacctcaaaagccagagatatagtccttagccgctattgcaattgatgcaactgaaattactgaaattagcactcctagtcctcatatttctcagtatagcaatgtaaatccatag