



Reference: PTXBC020614

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EX33; GPCR4; G-protein coupled receptor 84; inflammation-related G protein-coupled receptor EX33; G protein-coupled receptor 84		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 84
GENE ID	53831
ORIGIN SPECIES	Human
ACCESSIONS	BC020614
VECTOR	pDONR223
NCBI SYMBOL	GPR84
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1