

ZNF277-zinc finger protein 277 Gene

Reference: PTXBC020626

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NRIF4; ZNF277P; zinc finger protein 277; nuclear receptor-interacting factor 4; zinc finger protein (C2H2 type) 277; zinc finger protein 277 pseudogene		

For research use only.

GENE DESCRIPTION	zinc finger protein 277
GENE ID	11179
ORIGIN SPECIES	Human
ACCESSIONS	BC020626
VECTOR	pDONR223
NCBI SYMBOL	ZNF277
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggagaccgtgatggagctgcagcacagtcgggggtgtaggttatgggacagtaaggattgtatcctggagccgttccctgccagaaagtcagggtggcaccaccactttagaaggttccatctgtgcctgtattttgtggaagaacatttctgtggctgaacaagaca
aactctgaagcacatgattatgagcataagattgtcatagctgatgtcaagttggtgctgatttccaagggtacattttatattggaggaaagggttcactgaacagcccatcacagattttgtatgtaataagaattatccactgctccatttgaagaacaagaaattattttgtatgtga
cgttttaccagaagatagaattcttagagaagagcttcagaacagagactgagagaattctggaacaacagcagcaagaacgaatgataccaattttcatggcgtttgtatgtttgcaatgaagaattccttggaacagatctgtattttgaaccacatggccagagaacatgcttcaac
attggattgccagacaacattgtaactgcaatgaattttgtgtacattacagaaaagcttgacaatttgcagtcgtgtactgtgagaagaccttcagggacaaaaatacacttaagatcacatgaggaacacagcgcgttaagattatcctaagaacagagaatatgacagatttatg
tcataaattatttgtaaggctttttcataacttaaaatttgattgcagtactgcaaacaaaatcttctgaggaccaactaa