

C11orf1-chromosome 11 open reading frame 1 Gene

Reference: PTXBC020628

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0686 protein C11orf1; chromosome 11 open reading frame 1		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 1
GENE ID	64776
ORIGIN SPECIES	Human
ACCESSIONS	BC020628
VECTOR	pDONR223
NCBI SYMBOL	C11orf1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgctccagtggtctctgctcaaaattctctccagagacagaacctgcctgtttctcacaacccacactgtggcagcctgttaatgcagatggccatggtgaagtgtggacagattggaataatgtccaagtgtttccagtatggcgcgatgcaccactaatgagaata
cctattcaaacgtaccctgatggcgaactgaaccaggaaagatatgacctgagggaatatcgtagcagcccaaaccttgccttccagttggacactactttgaacaacatatgatacaagctacaacaacaaatgccactttcaacacatagatttaagcgagagcctcactggtccc
aggacatcaacctgaactggatcctccccgatacaaatgcacagaaaagtcaacttcatgaatagctattcaagccttaa