

ASAP1-ArfGAP with SH3 domain, ankyrin repeat and PH domain 1 Gene

Reference: PTXBC020631

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMAP1; CENTB4; DDEF1; PAG2; PAP; ZG14P; arf-GAP with SH3 domain, ANK repeat and PH domain-containing protein 1; 130 kDa phosphatidylinositol 4,5-bisphosphate-dependent ARF1 GTPase-activating protein; 130 kDa phosphatidylinositol 4,5-bisphosphate-dependent		

For research use only.

GENE DESCRIPTION	ArfGAP with SH3 domain, ankyrin repeat and PH domain 1
GENE ID	50807
ORIGIN SPECIES	Human
ACCESSIONS	BC020631
VECTOR	pDONR223
NCBI SYMBOL	ASAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgcacatctcagatgtgtgaagcatccattgcattttttattatttcttagttgttcttgacaaatttaaaacttttaaagattattcaagatga