



Reference: PTXBC020635

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FCNM; ficolin-1; M-ficolin; ficolin (collagen/fibrinogen domain-containing) 1; ficolin-A; ficolin-alpha; ficolin 1		

For research use only.

GENE DESCRIPTION	ficolin (collagen/fibrinogen domain containing) 1
GENE ID	2219
ORIGIN SPECIES	Human
ACCESSIONS	BC020635
VECTOR	pDONR223
NCBI SYMBOL	FCN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]