

RRP15-ribosomal RNA processing 15 homolog (S. cerevisiae) Gene

Reference: PTXBC020641

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RRP15-like protein; CGI-115; KIAA0507; ribosomal RNA-processing protein 15; ribosomal RNA processing 15 homolog		

For research use only.

GENE DESCRIPTION	ribosomal RNA processing 15 homolog (S. cerevisiae)
GENE ID	51018
ORIGIN SPECIES	Human
ACCESSIONS	BC020641
VECTOR	pDONR223
NCBI SYMBOL	RRP15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatgtaactggagccgtagcgtcggtgctggaagacgagccacagacactctgatagtgaaaggaagctgtgacggaaaggaccattttatctgatgacgcaatagaagctgacagtgagggtgatgctgagccctgtgacaagaaaatgaaatgatggagaa
tcaagtgtgggactaatatggcgtggcgagatgctatggcctaaagtcctcaacaagaaaactcctgaaagtaaacctactattctgtcaaaaataagaagctggaaaggaaaaagaaagttaaagcaagaagactagagaaaaataaacagcgtgataagggctggagtgaggaa
atgatgtgcagagtaaaagccagatgtttgccaagacaagagacagagagaatcttcagagaattgcaacaaggggtgtgtgcaattatattatgttcagaacatcaaaaagatgttgatgaaaaggttaaggaaagctggaaagttctatgagaagcgtgctaagttgataatcaactg
tttccaagaaagatttcacagtggtttgagagggatggatggaagtacaatgagactgttcagcaggaagaaacaaagccaaacagactgaagtgaatcagaagaagccaggttgacgatcctacgtgatgtttcatggtggagcatctatgaaagactgggacaagga
aagtgtggccagatgacagcagaccagaatctgcaagtgaactctgatacataa