

MRPL36-mitochondrial ribosomal protein L36 Gene

Reference: PTXBC020642

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRIP1; L36mt; MRP-L36; PRPL36; RPMJ; 39S ribosomal protein L36, mitochondrial; BRCA1-interacting protein 1; mitochondrial ribosomal protein L36		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L36
GENE ID	64979
ORIGIN SPECIES	Human
ACCESSIONS	BC020642
VECTOR	pDONR223
NCBI SYMBOL	MRPL36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaatctttataaggaaaatgtgaacctctgctctatctcagtcgtcacacgggtgaagcctcgagccctctccacatttctatttgatccattcgaggtgcagccccgtggctgtggaaccggggcagcagtcgcgtcaacttctctcaccggcctcctgccccatctgtgct
gcgctggggttcaaaaacaagactgtccttaagaagcgctgcaaggactgttacctgtgaagagggcggtcggtgtacgtctactgtaaaacccatccgaggcacaagcagagacagatgtag