

CD302-CD302 molecule Gene

Reference: PTXBC020646

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD302 molecule; CD302 antigen; BIMLEC; CLEC13A; DCL-1; DCL1; C-type lectin domain family 13, member A; DEC205-associated C-type lectin 1; type I transmembrane C-type lectin receptor DCL-1		

For research use only.

GENE DESCRIPTION	CD302 molecule
GENE ID	9936
ORIGIN SPECIES	Human
ACCESSIONS	BC020646
VECTOR	pDONR223
NCBI SYMBOL	CD302
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataagcatacataatgaagaagaaatgcttttatactggatacttgaaaaagcaatgaaaggccagatgatacctactagcatgttttatgacacagatgatcgagtttcaagtgtttgataattcaaatatgacattgataagtgacagaccaagatgatgaggatttagtt
gacacctgtgtttctgcacatcaagacaggtgaatgaaaaaggaattgtgaagtgttctctgtggaaggaaactatgcaaaacagctatccatacaaaaggaaatattatcagataaccacatttaatatcagcattggtgattgctagcacggttaatttgacagttttgggagcaat
catttggtctgtacaaaaacattctgattctgtttcaccacagtttttcaaccgaccccaatcaccttataatgaagactgtgttttgtagttggagaagaaatgaatatcctgttcaatttgactaa