

MS4A4A-membrane-spanning 4-domains, subfamily A, member 4 Gene

Reference: PTXBC020648

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4SPAN1; CD20-L1; CD20L1; HDCME31P; MS4A4; MS4A7; membrane-spanning 4-domains subfamily A member 4A; CD20 antigen-like 1; Fc epsilon receptor beta subunit homolog; four-span transmembrane protein 1; membrane-spanning 4-domains, subfamily A, member 4; membr		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 4
GENE ID	51338
ORIGIN SPECIES	Human
ACCESSIONS	BC020648
VECTOR	pDONR223
NCBI SYMBOL	MS4A4A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatcagacctacagcagacattgcaggcctgaagaaacacctttctgctgccatgacaacatgcaaggaaatggaacaggccatgccagggctggccctggtgtgccccagctgggaaacatggctgtcatattcacatctgtgaaaggattgcaagagaagttctgaag
ggagaacccaaagtcttggggtgtgagattctgactgacctgatgagccttagcatgggaataacaatgatgtgtatggcatctaactatggaagtaacctatttccgtgtatcgggtacacaatttgggggtcagtaattgtttattttaggagcctgtcaattgcagcaggaaatt
agaactacaaaaggcctgggtccgaggtagtctaggaatgaataatcaccagctctgactggctgcatcagggatcttaatacaacacatttagcttggcgtttttatcattccatcacccctactgttaactactatggcaactcaataatgtcatgggactatgccatctaatgggtctggatggc
atggtgctccttaagtgtgctggaattctgcatgtgtcctctctgctttggatgaaagtgtctgtgtaccctggtgggtgtgttaattctgccatcacattctacatggcagaacagcatctccacaccacttaatgaggtttga