

MRPL35-mitochondrial ribosomal protein L35 Gene

Reference: PTXBC020651

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L35mt; MRP-L35; 39S ribosomal protein L35, mitochondrial; mitochondrial ribosomal protein L35		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L35
GENE ID	51318
ORIGIN SPECIES	Human
ACCESSIONS	BC020651
VECTOR	pDONR223
NCBI SYMBOL	MRPL35
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgctctgcctttgctggtgcagtgagagcagcttcaggaatcctacggtccctgaatatttgacatctcaacacaccgcaactgtgtcaagaatgcctctcttatttctgcatgtccactggacgttttagtcataatcagacaccagttgttctccactcccagacttaccacatctga
gagaaacctgacatgtgggcataacctcagtgatccttaataagaatggccccgtgctccaagtgtcctgaagctgccagtcagatcttaacatactcagtgcaagaaaaggcaagagaaagaccgtgaaagctgtcatcgataggtttctcgactcattgtgaccttgggtgaggaga
aaggctggctataagaaaaattatgaaaaagacacctgcaaggagaagcgattgaggggaattgtattctgcaataaaaccagagtaaacctcttagataaaatgacgacgtcctctggaaggcgaaactga