

RAB8B-RAB8B, member RAS oncogene family Gene

Reference: PTXBC020654

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB8B, member RAS oncogene family; ras-related protein Rab-8B; RAB-8b protein		
For research use only.			

GENE DESCRIPTION	RAB8B, member RAS oncogene family
GENE ID	51762
ORIGIN SPECIES	Human
ACCESSIONS	BC020654
VECTOR	pDONR223
NCBI SYMBOL	RAB8B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaagacgtacgattatctctcaagctcctgctgacggcgactcggggtaggcaagacgtgcctcctgttccgcttctcagaggacgcttcaaccaccttccaccatcggaattgattttaaattagaacgatagaactagatggaagaaaattaagctcagatatg
ggacacagcgggtcaggaaagattccgaacaatcacgacgctactacagaggagccatgggcattatgctgctatgacatcacaaatgaaaaatcctttgacaattataaaattggatcagaacattgaagagcatgcctctccgatgtcgaagaatgatcctgggtaacaaatg
tgatatgaatgacaaaagacaagtgtcaaaagaaagaggggagaagctagcaattgactatgggattaaattcttgagacaagcgcaaaatccagtgcacaaatgtagaagaggcattttttacactgcacgagataatgacaaaactcaacagaaaaatgaatgacagcaattcagcag
gagcaggtggaccagtgaataaacagaaaaccgatcaagaagaccagttcttctgctcgtactttga