



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SPTLC3-serine palmitoyltransferase, long chain base subunit 3 Gene

Reference: PTXBC020656

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C20orf38; LCB2B; LCB3; SPT3; SPTLC2L; dJ718P11; dJ718P11.1; hLCB2b; serine palmitoyltransferase 3; LCB 3; SPT 3; long chain base biosynthesis protein 2b; long chain base biosynthesis protein 3; serine palmitoyltransferase, long chain base subunit 2-like (

For research use only.

GENE DESCRIPTION	serine palmitoyltransferase, long chain base subunit 3
GENE ID	55304
ORIGIN SPECIES	Human
ACCESSIONS	BC020656
VECTOR	pDONR223
NCBI SYMBOL	SPTLC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctaacccctggagggtggtgctgtttgcaacgggaaacttcaaacacaaagaacagagcaatggctcacaagaacgaaactgcacaaaagatggaatgttgaaaggagcccgacaaatgggaagccacattttatgataagctcattgttgaaatcgtttgagggaagcaccctctca
tgatggttttcaactcagggatgataaattggaacccctgtgttgcattctcagagactttttaaagaacatggggaataagaaatgcaacgcagctgttggaacgaaagaaacaaaagatcgtatgcgcaccctccgttgctgtttctcgaactcgtctctctaaagttgttccagaatggtg
tgatcaggtgctgagattcgaagaacatgagctgctcaggaaacttctcaggcgaacaaatcttcatctccactctgaataccccattaccgctctatggacgtgttttcaaacactctcctggcgcaactctatctcgcacatcaacaaactaa