

CXorf48-chromosome X open reading frame 48 Gene

Reference: PTXBC020662

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf48; BJHCC20A; cancer/testis antigen 55; BRCA2-interacting protein; tumor antigen BJ-HCC-20		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 48
GENE ID	54967
ORIGIN SPECIES	Human
ACCESSIONS	BC020662
VECTOR	pDONR223
NCBI SYMBOL	CXorf48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaggcttctgagactgtcttggccttctacggaggacggccgacctgcagagcgacagggccacagcagcagggcccccacaggtgacaccagttgacaactgtgaggaggtgtcacaagttctgtggtgattatggcatgattgatgagtcgatctactcagtagt
gatgtgtgactggcaacgtgcctctaaaagtggacaaaaagtaattggtgtggaagaagataaacacattatggattgagagcaatcaaggtggatgtgtgcctcgccatcttattgtgctggacctcagactcaggaaccagagtttaattggatgtgttactctataaatgaag
ataatatttatattagtaacagcatttattttccatagccattgtttcgaagattttgtgcctataaaggggtgactgttagaagttgaattatccactgagccagcattcacaacatcaaggcaacttctgtgaagccatccgtgtattcatagcgaagaggtctgcattactagcgtacatggaa
gaaatgggtgatagattatactattttttcaccttggattctgtgaaacttctgatggatagctacccaagtagatgatctgcaatgtggtcatggtggagagcattcagttctgttatttggagagcgatttctatcaccacgtgcataaatccagcagcggattccaggtacggag
gcctgggacggcccaaaagggaacgtcgagccaagcattta