

VPREB3-pre-B lymphocyte 3 Gene

Reference: PTXBC020666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	8HS20; N27C7-2; pre-B lymphocyte protein 3; pre-B lymphocyte gene 3; pre-B lymphocyte 3		

For research use only.

GENE DESCRIPTION	pre-B lymphocyte 3
GENE ID	29802
ORIGIN SPECIES	Human
ACCESSIONS	BC020666
VECTOR	pDONR223
NCBI SYMBOL	VPREB3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctgccggtgcctcagctcctctgatgggaccttctgtcagttccagacagtcctggccagctggatgcactgctgtcttccaggccaagtggctcaactctcctgcacgctcagcccccagcagtcaccatcagggactacgggtgtcctggtaccagcagcgggc
aggcagtgccctcgatatctcctactaccgctcggaggagatcaccaccggcctgtgacatccccgatcgattctcggcagccaaggatgaggcccaaatgcctgtgtcctcaccattagtcctgtgcagcctgaagacgacgcggattactactgctctgttgctacgcttta
gtccctag