

## MS4A7-membrane-spanning 4-domains, subfamily A, member 7 Gene

Reference: PTXBC020673

|                |                                                                                                                                                                                                                                                                |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | 4SPAN2; CD20L4; CFFM4; MS4A8; membrane-spanning 4-domains subfamily A member 7; CD20 antigen-like 4; CD20/Fc-epsilon-RI-beta family member 4; four-span transmembrane protein 2; high affinity immunoglobulin epsilon receptor beta subunit; membrane spanning |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | membrane-spanning 4-domains, subfamily A, member 7                                                |
| <b>GENE ID</b>          | 58475                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC020673                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | MS4A7                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgctattacaatcccaacatgggggtttctcacagctttacacaaaggcatcactatccctcaagagagaaacctggacacatgtacaaaacgaagattacctgcagaacgggctccaacagaaccaccgttctgggactgtccagatcctgtgttgccgttgattcaagct  
gggggccatctgtgtttgtccctacccctccactcaatccagcaattccaccacttggatgtctgggtacccatttttaggagctctgtgtttggcattactggatccctcattatctctggaacaaatcaactaagccctttgacctgagcagcttgacctcaatgcagtgcgtctgt  
tactgcaggagcaggcctctctcctgtgacagcatggtagccctgaggactgcctcacaattgtggctcagaatggattatctatcctcattgccttattcggagtactattccaatatgaaatcaagattgtctcctgaccagtgtcagtttaacagggtgcctagtgtgatgctc  
atcttcactgtctggagctcttattagctgcatacagttctgtctttggggaacagctctactccaacaacctgggagttcatttctcgaccagtcacaagatcatatccaacaggtcaaaaagagttcttcacggtcttgatataa