

## TNFSF13B-tumor necrosis factor (ligand) superfamily, member 13b Gene

Reference: PTXBC020674

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BAFF; BLYS; CD257; DTL; TALL-1; TALL1; THANK; TNFSF20; TNLG7A; ZTNF4; tumor necrosis factor ligand superfamily member 13B; ApoL related ligand TALL-1; B-cell-activating factor; B-lymphocyte stimulator; Delta4 BAFF; TNF and ApoL-related leukocyte expressed |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | tumor necrosis factor (ligand) superfamily, member 13b                                            |
| <b>GENE ID</b>          | 10673                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC020674                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TNFSF13B                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggatgactccacagaaaggagcagtcacgcttactcttgccttaagaaagagaagaatgaactgaaggagtgtgttccatctccacggaggaaagccctctgtccgatcctcaagacggaaagctgtgctgcaacctgtctgtgacctgtctgttgcctc  
acgggtgtgtcttaccaggtggccgctgcaaggggacgtggccagctccggcgagagctgcaggccaccacgcggagaagctgccagcaggagcagagccccaaggccgctggaggagctccagctgtcaccgcccgaactgaaaatcttgaaccacagctc  
caggagaaggcaactccagtcagaacagcagaataagcgtgcccgttcagggtccagaagaacagtcactcaagactgctgcaactgattgcagacagtgaaacacaaactatacaaaaaggatcttacacattgttccatggctctcagctttaaagggaagtgcctagaaga  
aaaagagaataaataattgtcaagaaactgttacttttatatattgtcaggtttatatactgataagacacacagcctgggacatctaattcagaggagaaggtccatgtcttggggatgaattgagtcgtgactttgttcgatgtattcaaatatgcctgaacactaccaataatt  
cctgctattcagctggcattgcaaaactggaagaggagatgaactcaacttgcataccaagagaaatgcacaaatcactggatggagatgtcacattttgtgctgattgaaactgctgtga