

HNMT-histamine N-methyltransferase Gene

Reference: PTXBC020677

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HNMT-S2; HNMT-S1; HMT; MRT51; histamine N-methyltransferase		

For research use only.

GENE DESCRIPTION	histamine N-methyltransferase
GENE ID	3176
ORIGIN SPECIES	Human
ACCESSIONS	BC020677
VECTOR	pDONR223
NCBI SYMBOL	HNMT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatcttccatgaggagctgttttctgaccacgggaatatgttgatcttccggaggtttctcaaccattccacgggaacaccagtgcatgcaggaaatcatggacaagaagctgccaggcataataggaaggattggagacacaaaatcagaaattaagattctagcatagggcgag
gtgcaggtgaaattgatcttcaaatctctccaaagttcaggctcaatccaggagtttctatcaacaatgaagttgtgagccaagtgctgaacaaattgccaatacaagagctgtgagccaagacatcgaacctcgagaacgtaaagtttgcctggcataaggagacatcatctgaatacc
aaagtagaatgttgagaaaaaggagcttcaaaagttgggactttatcatatgattcaaatgctgtattatgtaaaagacatccagctacccctgaaattctccatagctcttaggtaccaatgctaagatgctcattattgtgtgaggaagcagtggtgggacaagctgtgaaaaagtac
ggatcacgctttcccaggatgacctctgcccgtatcatcatcagatgacctcactcagatgctggacaacctagggccttaagtatgagtgctatgacctttgtccacatggatatactgactgctttattgatggtgatgaaaatggagacctgctttgggatttttactgaaacctgcaac
tttaatgccacagaccacctgatctcagagcagagcttgggaagatctacaagacctgaatttagtctgaagaaggagggaagggtcttttaataatactctgagttcatagtgattgagcgataa