



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

APOM-apolipoprotein M Gene

Reference: PTXBC020683

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

G3a; HSPC336; NG20; apo-M; NG20-like protein; alternative name: G3a, NG20; protein G3a

For research use only.

GENE DESCRIPTION	apolipoprotein M
GENE ID	55937
ORIGIN SPECIES	Human
ACCESSIONS	BC020683
VECTOR	pDONR223
NCBI SYMBOL	APOM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccaccaaatftgggcagctctctctacttctatggtattacttaacctactaccagtcgcttgagcagcagtcacacgaactctggcgctggatgggaaggaggtcccgagggtccacttggccagtggtactttatcgaggggcagctccccaagaggaggttggcaa
cttttgacctgttggaacaatgtctcaatatggctgctgctgctgctcccgatgcagctccacctctgctgctacatccgcgatgaagatgggctctgtgctcccgaaatggatctaccacctgactgaaggagacacagatctcagaactgaaggcccgctgatgaagactgagc
cttttccagctcatccccagtggtgaatcatgctgaatgagacagcgggtttaccagcgttttctctctacaatcgtccacacatctcccgaaagtgtgtggaggaattcaagtcctgactctgctgactcgaagcctcttattgactctaggaatcaagagccctgtgagc
tcttcccaatgca