



Reference: PTXBC020684

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	sushi repeat-containing protein SRPX; DRS; ETX1; HEL-S-83p; SRPX1; epididymis secretory sperm binding protein Li 83p; sushi-repeat-containing protein, X chromosome; sushi repeat containing protein, X-linked		

For research use only.

GENE DESCRIPTION	sushi-repeat-containing protein, X-linked
GENE ID	8406
ORIGIN SPECIES	Human
ACCESSIONS	BC020684
VECTOR	pDONR223
NCBI SYMBOL	SRPX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]