

LY96-lymphocyte antigen 96 Gene

Reference: PTXBC020690

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ESOP-1; MD-2; MD2; ly-96; lymphocyte antigen 96; myeloid differentiation protein-2; protein MD-2		

For research use only.

GENE DESCRIPTION	lymphocyte antigen 96
GENE ID	23643
ORIGIN SPECIES	Human
ACCESSIONS	BC020690
VECTOR	pDONR223
NCBI SYMBOL	LY96
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtaccattctgtttttccaccctgtttctcatalactgaagctcagaagcagfattgggtctgcaactcatccgatgcaagtattcatacacctactgtgataaaatgcaataccaattcaatgatgtaaccctgtatagaattgaaagatcacaaggattatgcacatttctaca
ttccaaggagagatttaaagcaattatattcaatctctataaactgtcaacacatgaatcttccaaagcgaagaagtattgcccaggatctgatgacgattacttttgcagagctctgaaggagagactgtgaatacaacaatcattctccccaagggaataaaatttctaagg
aaaatacaaatgtgttgaagctattctgggagccagaagaatgctctttgcttggagttgtcatcctacaccaacttaattcaaatag