

## DLEU1-deleted in lymphocytic leukemia 1 (non-protein coding) Gene

Reference: PTXBC020692

|                |                                                                                                                    |                     |                |
|----------------|--------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BCMS; DLB1; DLEU2; LEU1; LEU2; LINC00021; NCRNA00021; XTP6; deleted in lymphocytic leukemia 1 (non-protein coding) |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | deleted in lymphocytic leukemia 1 (non-protein coding)                                            |
| <b>GENE ID</b>          | 10301                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC020692                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | DLEU1                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagacctgtatattgatacagctgcaatttaaaaccgcctgccgctgttagagctttgccgttctccagcgtttacaggggttatcgacttaagcctcggacaactttaccagtgattctaccagaagggaatgaagaacagaaccttcaggaattgagtcacaatgcagacaatat  
caaatgggagattgttgcagggaagagattgatagtagtatttttactag