

CCL20-chemokine (C-C motif) ligand 20 Gene

Reference: PTXBC020698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CKb4; Exodus; LARC; MIP-3-alpha; MIP-3a; MIP3A; SCYA20; ST38; C-C motif chemokine 20; CC chemokine LARC; beta chemokine exodus-1; chemokine (C-C motif) ligand 20; liver and activation-regulated chemokine; macrophage inflammatory protein 3 alpha; small ind		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 20
GENE ID	6364
ORIGIN SPECIES	Human
ACCESSIONS	BC020698
VECTOR	pDONR223
NCBI SYMBOL	CCL20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgtaccaagattgtctcctgctgtcttgcagtcgtctctactccacctctgcggcgaatcagaagcaacttgcagtcgtctctggatcacagaccgtattcttcctctaaattattgtgggttcacacggcagctggccaatgaaggctgtgacatcaatgctatcatcttcacacaagaaaaagttgtctgtgtgcgcaaatccaaacagacttgggtgaatatattgtgcgtctcctcagtaaaaaagtcagaacatgtaa